

OUR CUSHIONS & POSITIONING SUPPORTS

SYST'AM® P314C / MIXTE CUSHION

MOULDED CUSHION MADE OF VISCOELASTIC GEL AND HR-FOAM

MATERIALS

40 kg/m³



HR foam



Visco Gel



Removable Polymaille® cover



Removable Polymaille® HD cover



NF EN ISO 597 - 1 & 2

Foam maintenance :



(Do not immerse in water)

Cleaning of the cover :



(UP TO 90°C)

Cleaning POLYMAILLE® HD cover :



< 1000 ppm

INDICATIONS

Risk

LOW




VERY HIGH

P314C / MIXTE-Cushion can be used, when :

- Seated position is necessary resulting from diseases or handicap and there is at the same time an increased risk of pressure sore.

SYST'AM® P314C (without cover)



SYST'AM® P314C (with cover)



NEW Polymaille® HD Cover (with welded seams)





Item code	Designation	Sizes (cm / inches) (W x L)	Patient weight (kg / lbs)	
			Min.	Max.
P314C40401HW	P314C - MIXTE® cushion with POLYMAILLE® cover	40 x 40 / 16 x 16	30 / 66	90 / 198
P314C45431HW			40 / 88	120 / 264
P314C40401HF	P314C - MIXTE® cushion with POLYMAILLE® HD cover	40 x 40 / 16 x 16	30 / 66	90 / 198
P314C45431HF			40 / 88	120 / 264

PRODUCT «ADVANTAGES»

COMBINATION OF GEL AND FOAM

- ↳ High-reticulated viscoelastic polyurethane gel : avoids gel migration thus preventing all punching of the support.
- ↳ The addition of a high-resilience foam base enables the cushion to adapt better to bony protuberances, thus providing a larger contact surface, therefore more effective distribution of pressures. Moreover, it improves patient comfort and stability.
- ↳ Because the two layers of the material are moulded together (≠ stuck), they are completely bonded, thus eliminating all risks of separating over time.



100% IMPERMEABLE

- ↳ Completely adapted to incontinent patients.



WARRANTY

2 YEARS

32

SYST'AM®

ZAE Route de Casseneuil - Le Lédat
47300 Villeneuve-sur-Lot - France

International Phone : +33 5 53 40 50 30
International Fax : +33 5 53 40 58 38

E-mail : info@system.com
www.system.com

