

DIAGNOSTIC[®]

PRO MESH

MESH NEBULIZER




**THE LEADING
RESPIRATORY
SPECIALIST**





INNOVATIVE MESH TECHNOLOGY



PARTICLE SIZE MMAD

3,8 μm

PORTABLE DEVICE (SMALL & LIGHT)

- possibility of using at any time, any place, even during the trip

ULTRA QUIET

- eliminates the fear of nebulization among younger children

UNIVERSAL

- for use with all medicines intended for nebulization

FULLY EQUIPPED

- universal model for the whole family: pediatric and adult masks, mouthpiece, power adapter

RECOMMENDED AT

cold

asthma

allergy

pulmonary
rehabilitation

bronchitis

COPD (Chronic Obstructive
Pulmonary Disease)

pneumonia



SPECIFICATIONS

Particle Size	MMAD 3,8 µm
Nebulization Rate	> 0,25 ml/min
Recommended fill volume	maximum 8 ml minimum 2 ml
Detects the absence of the drug in a container	LED orange LED indicator (POWER) turn-on for 8 sec, then power off
The amount of medicine that remains in the container after nebulization	< 0,3 ml
Power supply	3V DC ("AA" 1.5V alkaline battery x 2) AC adapter, Input: 100-240VAC, 0/60Hz, 0.2A Output: 5V DC, 1A
Power consumption	Approx. 1.2W
Vibrating Frequency	Approx. 113 kHz
Content package	Spray nozzle cap, mouthpiece, pediatric mask, adult mask, spray nozzle with medication chamber, alkaline batteries, case, instructions manual, power adapter



Manufacturer:
Diagnosis S.A.
ul. Gen. W. Andersa 38A
15-113 Białystok, POLAND
www.diagnosis.pl