



SYST'AM® P311C / GEL CUSHION

GEL CUSHION MADE OF VISCOELASTIC GEL

MATERIALS



Visco Gel



Removable POLYMAILLE® cover



Removable POLYMAILLE® HD cover



Cleaning of the gel:



(Do not immerse in water)

Cleaning POLYMAILLE® cover:




Cleaning POLYMAILLE® HD cover:



INDICATIONS

Risk





P311C / Gel cushion can be used, when:

- Seated position is necessary resulting from diseases or handicap and there is at the same time an increased risk of pressure sore.

VISCOELASTIC HIGH-RETICULATED POLYURETHANE GEL

→ Avoids gel migration and prevents the bottoming effect of the bony protuberances.



100% IMPERMEABLE
→ Completely adapted to incontinent patients.

AVAILABLE VERSIONS





SYST'AM® P311C /Gel Cushion
With POLYMAILLE® cover.



SYST'AM® P311C /Gel Cushion
With POLYMAILLE® HD cover
(with welded seams).



SYST'AM® P313C /Gel Cushion
Gel-Cushion with a removable
central pad.

DESIGNATION	ITEM CODE	 (W x L - cm / inches)	 min - MAX (kg / lbs)
SYST'AM® P311C with POLYMAILLE® cover	P311C40401HW	40 x 40 cm / 16 x 16"	40 - 120 kg / 88 - 264 lbs
	P311C45431HW	45 x 43 cm / 18 x 17"	40 - 140 kg / 88 - 308 lbs
SYST'AM® P311C with POLYMAILLE® HD cover	P311C40401HF	40 x 40 cm / 16 x 16"	40 - 120 kg / 88 - 264 lbs
	P311C45431HF	45 x 43 cm / 18 x 17"	40 - 140 kg / 88 - 308 lbs
SYST'AM® P313C with a removable central pad.	P313C40401HW	40 x 40 cm / 16 x 16"	40 - 100 kg / 88 - 220 lbs

