

# OUR CUSHIONS & POSITIONING SUPPORTS

## SYST'AM® P311C / GEL CUSHION

GEL CUSHION MADE OF VISCOELASTIC GEL

### MATERIALS



Visco Gel



Removable Polymaille® cover



Removable Polymaille® HD cover



NF EN ISO 597 - 1 & 2

Cleaning of the gel :



(Do not immerse in water)

Cleaning POLYMAILLE® cover :



(UP TO 90°C)

Cleaning POLYMAILLE® HD cover :



< 1000 ppm

### INDICATIONS

Risk



P311C can be used, when :

- Seated position is necessary resulting from diseases or handicap and there is at the same time an increased risk of pressure sore.

SYST'AM® P311C (without cover)



SYST'AM® P311C


↳ Gel-Cushion with a removable central pad.



SYST'AM® P311C (with cover)



**NEW** Polymaille® HD Cover (with welded seams)





Item code	Designation	Sizes (cm / inches) (W x L)	Patient weight (kg / lbs)	
			Min.	Max.
P311C40401HW	P311C - SYST'AM® gel with POLYMAILLE® cover	40 x 40 / 16 x 16	40 / 88	120 / 264
P311C45431HW				
P311C40401HF	P311C - SYST'AM® gel with POLYMAILLE® HD cover	40 x 40 / 16 x 16	40 / 88	120 / 264
P311C45431HF				
P313C40401HW	P313C - SYST'AM® gel with a removable central pad.	40 x 40 / 16 x 16	40 / 88	100 / 220

### PRODUCT «ADVANTAGES»

#### VISCOELASTIC HIGH-RETICULATED POLYURETHANE GEL

↳ Avoids gel migration and prevents the bottoming effect of the bony protubances.



#### 100% IMPERMEABLE

↳ Completely adapted to incontinent patients.

WARRANTY

2 YEARS

30

**SYST'AM®**

ZAE Route de Casseneuil - Le Lédât  
47300 Villeneuve-sur-Lot - France

International Phone : +33 5 53 40 50 30  
International Fax : +33 5 53 40 58 38

E-mail : info@system.com  
www.system.com

